VIC TO

Adjustable Artificial Sphincter

REPLACEMENT KITS

INFORMATION FOR SURGEONS - INSTRUCTIONS FOR USE

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DANSK ■

The logo for VICTO features the word "VICTO" in a bold, sans-serif font. Above the letter "I" is a stylized, grey, wing-like graphic element that points upwards and to the right.

VICTO

Adjustable Artificial Sphincter

GENERAL INFORMATION

INDICATIONS FOR USE

VICTO and VICTO+ are adjustable artificial urinary sphincters (AUS) used to treat male stress urinary incontinence caused by intrinsic sphincter deficiency (ISD) in cases such as post-prostatectomy incontinence.

The Replacement kits are intended to be used as replacement of cuff, pump or pressure regulating balloon of VICTO and VICTO+.

KIT CHARACTERISTICS

Cuff Replacement Kit (Ref.: KIT-RC-37-10-06; KIT-RC-40-10-06; KIT-RC-45-10-06; KIT-RC-50-10-06).

This kit contains the materials for replacement of the cuff of VICTO and VICTO+.

Balloon Replacement Kit (Ref.: KIT-RB-01)

This kit contains the materials for replacement of the Pressure Regulating Balloon of VICTO and VICTO+.

Pump Replacement Kit (Ref.: KIT-RP-03-10-06)

This kit contains the materials for replacement of the pump of VICTO and VICTO+.

DEVICE FUNCTION

The VICTO implant has been developed to simulate the normal function of the urethral sphincter and thereby maintain continence.

The device keeps the urethra closed by means of a urethral cuff, which prevents the passage of urine. The pressure in the cuff is regulated by a small silicone balloon (pressure-regulating balloon, or PRB). When the patient wants to void the cuff, it can be deflated by manual operation of the control pump, which is placed in the scrotum. The fluid in the cuff is transferred to the PRB, thus allowing the urine to pass through the urethra.

The cuff automatically reinflates via a slow-bleed valve in the control pump. Additionally, it can be manipulated for quick reinflation. The base of the control pump contains a self-sealing port, which enables the device to be interrogated percutaneously in order that extra fluid can be added to increase the pressure.

VICTO+ has an additional balloon placed in series with the cuff (stress-relief balloon, or SRB) to transfer additional fluid to the cuff during periods of high intra-abdominal pressure such as coughing and sneezing. This feature is called Conditional Occlusion.

The Replacement kits are used during VICTO/VICTO+ post-surgical replacement procedures of its parts, thereby allowing the surgeon to have flexibility once the product is implanted. This way, if replacement of the cuff/pump/balloon is necessary, the surgical procedure will be less invasive.

Components replacement might be due to an incorrect use/selection of the device, or product's component failure.

CONTRAINDICATIONS

The main contraindications for the implantation of the VICTO+ and VICTO are listed below:

1. Unmanaged detrusor instability;
2. Problems affecting manual dexterity or motivation, which may prevent the patient from operating the device;
3. Acute urinary tract infection which may lead to post-operative complications;
4. Patients whom the surgeon determines to be not suitable due to risks associated with open surgical procedures and/ or with the patients' medical history (physical or mental problems);
5. Known sensitivity to silicone rubber.

WARNINGS

Improper cuff sizing may limit the efficacy of therapy or cause urethral erosion.

Improper filling of components may limit the performance of the product or cause permanent damage.

All components are for single use. Thus, they should not be reused or resterilized, as this can potentially result in compromised device performance and increased risk of inappropriate resterilization and cross contamination.

The product must not be used if the package is open or damaged.

The product must not be used after the expiration date specified in the packaging.

PRECAUTIONS

Improper cuff sizing may result in tissue erosion, migration of components or continued incontinence.

Unsuccessful outcomes may result from improper sizing and/or filling of components.

The kits' components must be handled only by surgeons and/or operating-room staff who are familiar with the device and surgical technique.

Adequate precautions must be taken to avoid contamination during the intraoperative procedure.

Operating-rooms conditions should be in accordance with hospital, administrative and/or local government policy.

After use, the non-implanted product and packaging should be disposed in accordance with hospital, administrative and/or local government policy.

Right and complete connection between the connector and the tubing (perfect fit) should be checked. Confirm that there is no loss of fluid in the joints. Verify that the tubing is not kinked, either pre- or post-connector.

MAGNETIC RESONANCE IMAGING (MRI) ENVIRONMENT

Non-clinical testing performed in the worst-case scenario has demonstrated that VICTO/VICTO+ implant is MR conditional. A patient

with this device can be safely scanned in an MR system that meets the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, **with**
- **Maximum spatial field gradient of 12,800 G/cm (128 T/m)**
- Maximum force product of 231 T² /m
- **Theoretically estimated** maximum **whole-body** averaged (WBA) specific absorption rate (SAR) of
< 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the VICTO+/VICTO implant is expected to produce a maximum temperature rise of less than

1.0°C (**2 W/kg, 1.5 Tesla**) **RF-related temperature increase with a background temperature increase of ≈ 1.0°C (2 W/kg, 1.5 Tesla)**

0.9°C (**2 W/kg, 3 Tesla**) **RF-related temperature increase with a background temperature increase of ≈ 1.0°C (2 W/kg, 3 Tesla)**

after 15 minutes of continuous scanning.

In non-clinical testing, the image caused by the device extends approximately 18.3 mm from the implant at the port of the pump when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

POSSIBLE COMPLICATIONS

The use of VICTO/VICTO+ may result in complications associated with the methods used during the surgical procedure or reaction to the foreign body.

The efficacy of the device may be reduced in the patient who has been treated with radiotherapy for local recurrence of prostate cancer.

Some complications may require complete or partial removal of the device.

Infections that do not respond to antibiotic treatment and/or chronic pain related to the device failure may require the removal of the device.

Complications reported with this or other surgical procedures for the treatment of incontinence include:

- Urinary retention due to the patient's inability to empty the cuff correctly
- Fluid loss due to disruption of a connection or a pinhole leak
- Cuff and/or pump erosion (characterized by pain, hematuria, dysuria and recurrent incontinence)
- Abnormal detrusor function
- Urethral atrophy
- Bladder-neck or urethral closure
- Infection
- Mechanical failure
- Explantation
- Pain/discomfort
- Device malfunction
- Delayed wound healing
- Migration
- Recurrent incontinence
- Bladder spasm
- Difficult activation
- Tissue erosion
- Swelling
- Fistula formation
- Hematoma
- Hydrocele
- Tissue erosion/infection
- Patient dissatisfaction
- Positional incontinence
- Wound infection
- Urinary retention
- Difficult deactivation

Clinical literature reports that a small percentage of patients may continue to experience incontinence after an Artificial Urinary Sphincter placement.

COMPONENTS

PACKAGING

The Replacement Kits are supplied sterile and packaged in double plastic pouches to facilitate visualization of the device. Each kit is delivered in a labeled cardboard box for additional protection.

INDIVIDUAL COMPONENTS

THE CUFF REPLACEMENT KIT CONTAINS:

- 1 Replacement Cuff 3.7cm (RC-37), 4.0cm (RC-40), 4.5cm (RC-45), or 5.0 cm (RC-50)
- 2 connectors (STC-12)
- 1 blunt injection needle (BN-15G)
- 2 non-coring injection needles (NCN-25G)
- 1 syringe 6 ml (SYR-06)
- 1 syringe 10 ml (SYR-10)

THE BALLOON REPLACEMENT KIT CONTAINS:

- 1 Replacement Balloon (RB-01)
- 2 connectors (STC-12)
- 2 blunt injection needles (BN-15G)
- 2 non-coring injection needles (NCN-25G)
- 1 syringe 6 ml (SYR-06)
- 1 syringe 10 ml (SYR-10)

THE PUMP REPLACEMENT KIT CONTAINS:

- 1 Replacement Pump (RP-03)
- 3 connectors (STC-12)
- 2 blunt injection needles (BN-15G)
- 2 non-coring injection needles (NCN-25G)
- 1 syringe 6 ml (SYR-06)
- 1 syringe 10 ml (SYR-10)
- 1 protection tube (PT-01)

STERILIZATION – SUPPLIED STERILE

All components of the Replacement Kits are supplied sterile.

Sterility is maintained only if the packaging is intact and undamaged. Therefore, do not use the product if the package is open or damaged.

SHELF LIFE

For shelf life, refer to the label on the packaging of each individual components or the cardboard box label. Do not use after the expiration date specified in the packaging.

SINGLE USE

All components are intended for SINGLE USE ONLY. DO NOT REUSE or RESTERILIZE them, as this can potentially result in compromised device performance and increased risk of inappropriate re-sterilization and cross-contamination.

CARE AND STORAGE OF STERILIZED COMPONENTS

In order to protect the packaging integrity and sterile condition of the components, it is important that, during storage, they remain in the cardboard box on a protected shelf or inside a cabinet. The environment should be clean, dry and near room-temperature.

IMPLANTATION

PRE-OPERATIVE

SURGICAL TEAM PREPARATION

The surgeon and operating-room staff should be familiar with the operating procedure.

The members of the surgical team should scrub their hands for ten minutes using a povidone-iodine solution (provided the patient is not iodine sensitive) or another approved hospital scrub procedure.

PATIENT PREPARATION

The patient should be fully informed of the procedure and must give consent. The surgeon should discuss the possible risks associated with the surgery, anaesthetic and device.

Prior to surgery, prophylactic antibiotic therapy should be administered according to the hospital's approved procedures.

Once in the operating room, the patient is placed in the lithotomy position. Based on the component to be replaced, the abdominal and genital areas are shaved and cleansed with povidone-iodine solution or using an approved hospital preoperative scrub procedure.

The patient is prepped and draped for perineal and/or lower abdominal incision.

EQUIPMENT NEEDED

The implantation requires a standard urethroplasty surgical kit, the VICTO Accessories Kit (Ref.: KIT-AC-01), and various other supplies. Each surgeon may also have preferences about what equipment he/she requires.

Additional equipment is listed below:

- 12 Fr Foley catheter
- Shod forceps (e.g., hemostats)
- Sterile tray
- Sterile 50 ml syringe
- Sterile saline solution
- Umbilical tape

UNPACKING THE COMPONENTS

The integrity of both packages should be inspected. If there is a suspected breach of sterility the device should be returned to the distributor or manufacturer.

Select the replacement cuff after the correct measurement of the urethral size using the measurement tape provided in the Accessories Kit.

A number of labels detailing the batch number and use by date are supplied. These should be retained and placed in the appropriate

surgical and patient notes as established in the hospital procedure.

EQUIPMENT PREPARATION

PREPARING HEMOSTATS

During the procedure, it is useful to clamp some of the device tubes. In order to protect the device tubing from damage when clamping, cover the jaws of the hemostat with small sections of silicone tubing provided in Accessories Kit. All the teeth of the hemostat should be covered. The jaws should only be clamped to the first click to prevent excessive pressure on the tubing.

COMPONENT SET UP

CUFF

Before opening the pouches that keep the cuff sterile, perform the measurement of the urethra to ensure that the appropriate cuff size is selected.

Use of scissors to trim the cuff tubing length squarely to fit patient's anatomy.

Attach the blunt injection needle to the 10ml syringe.

Fill the syringe with approximately 5ml of sterile saline solution.

Hold the cuff with one hand and squeeze it until it is deflated.

Insert needle into the end of the cuff tubing.

Aspirate any air left in the cuff.

Fill the cuff with 1 ml to 3 ml, depending on cuff size.

NOTE: Do not over fill the cuff. Excessive fluid can stretch cuff material.

Squeeze the cuff and let it refill. Repeat it until most air bubbles are out of the cuff. Clamp off the cuff tubing 5cm below the needle.

Place the cuff into a sterile tray containing sterile saline solution.

BALLOON SET UP

Use scissors to trim the balloon tubing length squarely to fit patient's anatomy.

Attach the blunt injection needle to the 10ml syringe.

Fill the syringe with 10ml of sterile saline solution.

Hold the balloon with one hand and squeeze it until it is deflated.

Insert the needle into the end of the balloon tubing.

Fill the balloon with 3ml of sterile saline solution.

Aspirate any air left in the balloon.

Fill the balloon with 6ml of sterile saline solution.

NOTE: Do not overfill the balloon. Excessive fluid can stretch the balloon material.

Squeeze the balloon and let it refill. Repeat it until most of the air bubbles are out of the balloon. Clamp off the balloon tubing 5cm below the needle.

Place the balloon into a sterile tray containing sterile saline solution.

PUMP SET UP

Place end of each tubing in a basin with sterile saline solution. Use scissors to trim the control pump tubing length squarely to fit patient's anatomy.

Hold the control pump with the tubes on top.

With the tubing submerged, squeeze and release the control pump until all air in the control pump and tubing has been displaced with sterile saline solution.

In this position, clamp each tube 5cm from end using hemostats.

Place the control pump into a sterile tray containing sterile saline solution.

SURGICAL PROCEDURES

CUFF REPLACEMENT

Place an indwelling catheter (Foley 12) into the urethra to facilitate identification during the dissection.

Make a perineal incision and dissect around the VICTO+/VICTO cuff.

Place some umbilical tape around the urethra. Detach the cuff and released it to the urethra.

Place the measurement tape, supplied with the Accessories Kit, around the urethra where the replacement cuff is to be implanted. It should fit snugly without constricting the urethra. NOTE: Before measuring the urethra, take out the indwelling catheter.

Select the Cuff Replacement Kit that corresponds to the urethra size.

Using a hemostat, clamp the tubing of the cuff to be replaced, 5cm above the cuff tubing connector

Use scissors to trim the cuff tubing length squarely

Flush both pieces of tubing with sterile solution

Hold the titanium connector from the center and attach each tubing until achieve the central reference of the connector (as indicated in Figure 1).


FIGURE 1

Detach the hemostat and check the right and complete connection between the connector and tubing (perfect fit). Confirm that there is no loss of fluid in the joint. Verify that there is no kink in the tubing, either pre- or post-connector.

Squeeze the pump to empty the cuff completely.

The cuff is placed around the urethra and the tab is passed through the loop. Using a hemostat, the button is passed through the hole to secure the cuff in position. The cuff should not be tight (i.e., it should be able to rotate freely around the urethra).

The perineal incisions can be closed and dressed.

PRESSURE-REGULATING BALLOON REPLACEMENT

Make a lower abdominal incision to find the PRB.

Using a hemostat, clamp the tubing of the PRB to be replaced, 5cm above the PRB connector

Use scissors to trim the PRB tubing length squarely

Flush both pieces of tubing with sterile solution.

Hold the titanium connector from the center and attach each

piece of tubing until the central reference of the connector is achieved (as indicated in Figure 1).

Detach the hemostat and check for the proper connection between the connector and the tubing (perfect fit). Confirm that there is no loss of fluid in the joint. Verify that there is no kink in the tubing, either pre- or post-connector.

The abdominal incisions can be closed and dressed.

CONTROL PUMP REPLACEMENT

Make a lower abdominal incision to access the pump to be replaced.

Using two hemostat, clamp both pieces of tubing of the control pump, 5cm above the top of it.

Use scissors to trim the pump tubing length squarely.

Flush both pieces of tubing with sterile solution.

Hold the titanium connector from the center and attach each tubing until you achieve the central reference of the connector (as indicated in Figure 1)

Detach the hemostat and confirm a proper connection between the connector and tubing (perfect fit). Confirm that there is no loss of fluid in the joint. Verify that there is no kink in the tubing, either pre- or post-connector.

The abdominal incisions can be closed and dressed.

POST-OPERATIVE CARE

IMMEDIATELY POST-OPERATION

An antibiotic cover should be prescribed at the discretion of the surgeon.

For the situation of cuff replacement, an indwelling catheter should remain in place for 24 hours.

After 24 hours the indwelling catheter should be removed and the patient should try to void urgently.

If the patient is unable to void without a considerable amount of residual urine, it may be necessary to teach the patient intermittent catheterization in case of retention.

The patient should continue to use pads or condom drainage to contain any leakage.

In the case of pump or balloon replacement, a pressurization procedure must be made day after the replacement process.

DISCHARGE

The patient can be discharged 2 to 4 days postoperatively at the discretion of the surgeon and depending on his general health.

Antibiotic cover should be prescribed at the discretion of the surgeon.

An outpatient appointment should be made within 2 to 4 weeks following the surgical procedure.

In the case of cuff replacement, an appointment for pressurization should be made for 4 to 6 weeks following the cuff change.

PRESSURIZATION PROCEDURE

The VICTO+/VICTO implant is pressurized by injecting sterile saline solution through the self-sealing port in the base of the control pump. The pressure in the device depends upon the amount of fluid injected. This feature

allows the device to be set at any desired pressure without having to change the PRB.

EQUIPMENT REQUIRED

The pressurization procedure is best carried out as a sterile procedure in a clean clinical space.

- Sterile saline solution for injection
- Victo adjustment kit
 - 6ml syringe
 - Non-coring needle 25Gx5/8"
- Sterile dressing pack
- Sterile drapes
- Sterile gloves

IMPORTANT

It is vital that the wall of the pump unit is not pierced by the needle. If this occurs, the device may be permanently damaged, necessitating explantation of the entire device.

It is essential that only the 25G x 5/8" needle provided in the Adjustment Kit be used and that the needle is always aligned parallel to the sides of the pump so that it will penetrate the self-sealing port straight and through the "guard ring." Only the base of the pump with its metal guard ring may be injected safely.

PATIENT PREPARATION

- Explain the procedure to the patient.
- Lay the patient supine and identify the location of the control pump inside the scrotum.
- Cleanse the scrotum in accordance with the hospital's approved procedure.

NOTE: Local anesthesia around the area may be applied in accordance with the hospital's approved procedure.

PRESSURIZATION PROCESS

- Fill the syringe with a sterile saline solution for injection.
- Attach the non-coring needle 25G x 5/8" to the syringe.
- With syringe pointed upward, flush 1.0ml of sterile saline solution and expel excess air.
- Fix the self-sealing port of the control pump, using two fingers to pull the skin taut over the end of the control pump.
- Push the needle slowly through the skin into the center of the control pump self-sealing port. Just prior to reaching the control pump, clear the end of the needle by injecting 0.1 ml of sterile saline solution into the tissue space.
- Push the needle into the self-sealing port, ensuring that the needle remains parallel with the walls of the control pump.
- Activate the product, injecting 2ml of sterile saline solution.
- Remove the needle slowly.
- Instruct the patient to pump the device about 3 times to redistribute fluid. Rapid reinflation can be achieved by squeezing the top of the control pump.
- Perform the test of efficacy (Please refer to "Post-pressurization continence tests.")

POST-PRESSURIZATION CONTINENCE TESTS

TESTS OF EFFICACY

These can be done immediately following pressurization process. After showing the patient how to deflate the cuff for voiding, he can be asked to fill his bladder by drinking until a desire to void is reached.

The patient should be asked to perform 5 strong coughs, 5 seconds of Valsalva manoeuvre and, if possible, 5 squats. The volume leaked can be measured and recorded. The patient should then be asked to void into a flow machine in order that the flow rate and volume can be measured and recorded.

LONGER-TERM TESTS OF EFFICACY

The patient can be asked to continue wearing pads (or Conveen leg bag) to record pad use and/or leaked volume by weight.

PATIENT EDUCATION

Subsequent to device pressurization, the surgeon should show the patient how to pump the device in order to allow urination.

The rapid reinflation capability can also be demonstrated. If the patient is using intermittent catheterization to empty the bladder, the importance of pumping before catheter insertion should be stressed.

A flow-rate test should be performed to ensure that the patient is pumping down sufficiently.

An emergency card should be issued to inform healthcare professionals that the patient has an implant and that catheterization should only be performed after pumping down the device to deflate the cuff.

IMAGING

X-RAY VISIBILITY OF THE VICTO IMPLANT

The VICTO+/VICTO does not contain radio-opaque solution, and consequently it cannot be imaged using plain radiography (x-ray). However, there are 2 other methods that can be used to image the device in order to investigate function.

MAGNETIC RESONANCE

The VICTO/VICTO+ implant can be observed imaged using MRI (Magnetic Resonance Imaging) as described in the paper by Deng et al.

FLEXIBLE CYSTOMETRY

The function of the cuff can be visualized using flexible cystometry with the cuff inflated and deflated.

ULTRASOUND DIAGNOSTICS

The balloons of the VICTO/VICTO+ implant can be easily visualized using ultrasound. With appropriate settings, the outline, particularly the highlighted top and bottom of the balloons, can be seen.

The straight line distance between these highlights will give the diameter of the balloons. When the cuff is normally inflated between voidings, the SRB and PRB should be very nearly of the same diameter and pressure.

From the diameter of one of the balloons, and using the following formula, the added volume can be determined.

$$[(4/3\pi(D/2)^3)-9]$$

With this ultrasound technique it is also possible to determine the normal functioning of the device during voiding. Pump the sphincter cuff down with about 2 or 3 good squeezes and notice the enlargement of the PRB while the SRB diminishes in size. After waiting for approximately 10 to 15 minutes, during which time the cuff will be fully reinflated, the two balloons should have returned to approximately the same size.

BIBLIOGRAPHY & REFERENCES:

- Knight SL, Susser J, Greenwell T, Mundy AR, Craggs MD. A New Artificial Urinary Sphincter with Conditional Occlusion for Stress Urinary Incontinence: Preliminary Clinical Results. *European Urology*.2006; 50: 574–580.
- Deng J, Hall-Craggs MA, Craggs MD, Rickards D, Knight SL, Linney AD, Mundy AR. Three-dimensional MRI of the male urethrae with implanted artificial sphincters: initial results. *The British Journal of Radiology*. 2006;79.
- Craggs M, Knight S, Susser J, Kirkham A, Rickards D. Ultrasound imaging as a diagnostic tool for assessing a new artificial urinary sphincter. *Neurourology & Urodynamics*. Proceedings of the International Continence Society, 2006. 174.
- García Montes F, Knight SL, Greenwell T, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007; 31(7): 752-758.
- Hübner W, Oliver S, Schlarp O. Evaluation of the Flowsecure® Artificial Urinary Sphincter (AUS) Artificial Sphincter: Preliminary results. Proceedings of the International Continence Society, Rotterdam. 2007, 499.
- Hübner W. New developments in the area of hydraulic sphincters. *J.Urol.Urogynäkol*. 2007; 14:23-24
- García-Montes F, Vicens-Vicens A, Ozonas-Moragues M, Oleza- Simo J, Mora-Salvá A. Surgical implantation of the new Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter in the female bladder neck. *Urol Int*. 2007;79(2):105-10.
- García Montes F, Vicens Vicens A, Ozonas Moragues M, Pizá Reus P, Mora Salvá A, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007 Sep;31(8):872-879. [In Spanish].
- Khatib-Shahidi, K., Grein M. - Flowsecure® Artificial Urinary Sphincter (AUS) First experiences. Proceedings of NRW Urology, Congress,2008. HELIOS Hospital Schwelm, Germany Poster P60.
- Bestard Vallejo JE, García Montes F, Cecchini Rosell L, Mir Maresma C, Morote Robles J. Implantation technique of the artificial urinary sphincter "Flowsecure® Artificial Urinary Sphincter (AUS)" in the bulbar urethra. *Arch Esp Urol*.2009;62(3):195-200.
- Alonso Rodríguez D, Fernández Barranco L, García-Montes F, Mora Salvá A, Ozonas Moragues M. Follow-up of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter. *ActasUrol Esp*. 2009 Oct;33(9):952-5. Review.
- García-Montes F. Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter for the treatment of stress urinary incontinence after radical prostatectomy. *Arch Esp Urol*. 2009; 62(10): 845-850. [Spanish].
- Fernandez Barranco L, D. Rodríguez A, Garcia Montes F, Mora Salvá A, Ozonas Moragues M. MRI of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter following radical prostatectomy. *Arch. Esp. Urol*.2010; 63 (4): 313-314.
- Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. One hundred Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincters. *European Urology Suppl*. 2011; 10(2):309.
- Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. Four years experience with the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter: Problems and solution. *Neurourology & Urodynamics*. 2011 Vol Abstract 250.

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Knight SL, Susser J, Greenwell T, Mundy AR, Craggs MD. A New Artificial Urinary Sphincter with Conditional Occlusion for Stress Urinary Incontinence: Preliminary Clinical Results. *European Urology*.2006; 50: 574–580.

Deng J, Hall-Craggs MA, Craggs MD, Rickards D, Knight SL, Linney AD, Mundy AR. Three-dimensional MRI of the male urethrae with implanted artificial sphincters: initial results. *The British Journal of Radiology*. 2006;79.

Craggs M, Knight S, Susser J, Kirkham A, Rickards D. Ultrasound imaging as a diagnostic tool for assessing a new artificial urinary sphincter. *Neurourology & Urodynamics*. Proceedings of the International Continence Society, 2006. 174.

García Montes F, Knight SL, Greenwell T, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007; 31(7): 752-758.

Hübner W, Oliver S, Schlarp O. Evaluation of the Flowsecure® Artificial Urinary Sphincter (AUS) Artificial Sphincter: Preliminary results. Proceedings of the International Continence Society, Rotterdam. 2007, 499.

Hübner W. New developments in the area of hydraulic sphincters. *J.Urol.Urogynäkol*. 2007; 14:23-24

García-Montes F, Vicens-Vicens A, Ozonas-Moragues M, Oleza- Simo J, Mora-Salvá A. Surgical implantation of the new Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter in the female bladder neck. *Urol Int*. 2007;79(2):105-10.

García Montes F, Vicens Vicens A, Ozonas Moragues M, Pizá Reus P, Mora Salvá A, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007 Sep;31(8):872-879. [In Spanish].

Khatib-Shahidi, K., Grein M. - Flowsecure® Artificial Urinary Sphincter (AUS) First experiences. Proceedings of NRW Urology, Congress,2008. HELIOS Hospital Schwelm, Germany Poster P60.

Bestard Vallejo JE, García Montes F, Cecchini Rosell L, Mir Maresma C, Morote Robles J. Implantation technique of the artificial urinary sphincter "Flowsecure® Artificial Urinary Sphincter (AUS)" in the bulbar urethra. *Arch Esp Urol*.2009;62(3):195-200.

Alonso Rodríguez D, Fernández Barranco L, García-Montes F, Mora Salvá A, Ozonas Moragues M. Follow-up of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter. *Actas Urol Esp*. 2009 Oct;33(9):952-5. Review.

García-Montes F. Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter for the treatment of stress urinary incontinence after radical prostatectomy. *Arch Esp Urol*. 2009; 62(10): 845-850. [Spanish].

Fernandez Barranco L, D. Rodríguez A, Garcia Montes F, Mora Salvá A, Ozonas Moragues M. MRI of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter following radical prostatectomy. *Arch. Esp. Urol*.2010; 63 (4): 313-314.

Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. One hundred Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincters. *European Urology Suppl*. 2011; 10(2):309.

Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. Four years experience with the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter: Problems and solution. *Neurourology & Urodynamics*. 2011 Vol Abstract 250.

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Sphinx Medical Limited,