ARTIFICIAL URINARY SPHINCTER (AUS)

INFORMATION FOR SURGEONS - INSTRUCTIONS FOR USE

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DANSK ■

The logo for VICTO features the word "VICTO" in a bold, sans-serif font. Above the letter "I" is a stylized, grey, wing-like graphic element that points upwards and to the right.

VICTO

Adjustable Artificial Sphincter

GENERAL INFORMATION

INDICATIONS FOR USE

VICTO is an adjustable artificial urinary sphincter (AUS) used to treat male stress urinary incontinence caused by intrinsic sphincter deficiency (ISD) in cases such as post-prostatectomy incontinence.

DEVICE CHARACTERISTICS

STARTER KIT (KIT-AUS-37-WS; KIT-AUS-40-WS; KIT-AUS-45-WS; KIT-AUS-50-WS)

The VICTO implant is a medical grade silicone device which consists of a pressure regulating balloon (PRB), an urethral cuff, and a control pump with a self-sealing port.

The necessary materials to set up the device before to be implanted are included in the kit.

ACCESSORIES CHARACTERISTICS

ACCESSORIES KIT (KIT-AC-01)

The Accessories Kit contains materials that are necessary for the implantation procedure.

ADJUSTMENT KIT (KIT-AD-06)

This kit contains the necessary materials to adjust the system internal pressure by adding or extracting sterile saline solution.

DEVICE FUNCTION

The VICTO implant has been developed to simulate the normal function of the urethral sphincter in order to maintain continence.

The device keeps the urethra closed by means of an urethral cuff which prevents the passage of urine. The pressure in the cuff is regulated by a small silicone balloon

(pressure regulating balloon - PRB). When the patient wishes to void the cuff, it can be deflated by manual operation of the control pump which is placed in the scrotum. The fluid in the cuff is transferred to the PRB allowing the urine to pass through the urethra.

The cuff automatically re-inflates through a slow bleed valve in the control pump; in addition, it can be manipulated for fast re-inflation. The base of the control pump contains a selfsealing port which enables the device to be interrogated percutaneously in order that extra fluid can be added to increase the pressure.


Figure 1


Figure 2

A) The cuff is deflated by squeezing the control pump unit located in the scrotum; fluid is transferred to the PRB, leaving the cuff empty and allowing the patient to void.

B) The cuff slowly re-inflates through a slow bleed valve. Rapid re-inflation can be achieved by squeezing the top of the control pump. The cuff exerts sufficient pressure to maintain continence without compromising the urethral tissues.

CONTRAINDICATIONS

The main contraindications for the implantation of VICTO are listed as follows:

1. Unmanaged detrusor instability.
2. Problems affecting manual dexterity or motivation, which may prevent the patient from operating the device.
3. Acute urinary tract infection which may lead to post-operative complications.
4. Patients whom the surgeon determines to be not suitable due to risks associated with open surgical procedures and/or with the patients' medical history (physical or mental problems).
5. Known sensitivity to silicone rubber.

WARNINGS

Improper cuff sizing may cause low efficacy therapy or urethral erosion.

Improper filling of components may cause low performance of the product or implant permanent damage.

All components are for single use. Therefore, they should not be re-used or re-sterilized, as this can potentially result in compromised device performance and increased risk of inappropriate re-sterilization and cross contamination.

The product must not be used if the package is open or damaged.

The product must not be used after the expiry date specified on the packaging.

PRECAUTIONS

Improper cuff sizing may result in tissue

erosion, migration of components, or continued incontinence.

Unsuccessful outcomes may result from improper sizing and/or filling of components.

The kits' components must be handled only by surgeons and/or operating room staff who are familiar with the device and the surgical technique.

Adequate precautions must be taken to avoid contamination during the intraoperative procedure.

Operating rooms conditions should be in accordance with hospital, administrative and/or local government policy.

After use, the non-implanted product and packaging should be disposed in accordance with hospital, administrative and/or local government policy.

MAGNETIC RESONANCE (MR) ENVIRONMENT

Non-clinical testing performed in the worst case product has demonstrated that VICTO implant is MR conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla, with
- Maximum spatial field gradient of 12,800 G/cm (128 T/m)
- Maximum force product of 231 T²/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of
< 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the VICTO implant is expected to produce a maximum temperature rise of less than 1.0°C (2 W/kg, 1.5 Tesla) RF-related temperature increase with a background temperature increase of ≈ 1.0°C (2 W/kg, 1.5 Tesla) 0.9°C (2 W/kg, 3 Tesla) RF-related temperature increase with a background temperature increase of ≈ 1.0°C (2 W/kg, 3 Tesla)

after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 18.3 mm from the VICTO implant at the “port of the pump” when imaged with a gradient echo pulse sequence and a 3 Tesla MR system.

POSSIBLE COMPLICATIONS

The use of VICTO may result in complications associated with the methods used during the surgical procedure or reaction to the foreign body.

The efficacy of the device may be reduced in patients who had been treated with radiotherapy for local recurrence of prostate cancer.

Some complications may require complete or partial removal of the device.

Infections that do not respond to antibiotic treatment, chronic pain related with the device failure may require the removal of the device.

Complications reported with this or other surgical procedures for the treatment of incontinence include:

- Urinary retention due to inability of the patient to empty the cuff correctly
- Fluid loss due to disruption of a connection or a pinhole leak
- Cuff and/or pump erosion (characterized by pain, haematuria, dysuria, and recurrent incontinence)
- Abnormal detrusor function
- Urethral atrophy
- Bladder neck / urethral closure
- Infection
- Mechanical failure
- Explantation
- Pain/discomfort
- Device malfunction
- Delayed wound healing
- Migration
- Recurrent incontinence
- Bladder spasm
- Difficult activation
- Tissue erosion

- Swelling
- Fistula formation
- Haematoma
- Hydrocele
- Tissue erosion/infection
- Patient dissatisfaction
- Positional incontinence
- Wound infection
- Urinary retention
- Difficult deactivation

Clinical literature reports that a small percentage of patients may continue experiencing incontinence after an Artificial Urinary Sphincter placement.

COMPONENTS

PACKAGING

The Starter kit, Accessories kit, and Adjustment kit are supplied sterile and packaged in double plastic pouches to enable visualization of the device. Each kit is delivered in a labeled cardboard box for additional protection.

INDIVIDUAL COMPONENTS

STARTER KIT

The VICTO implant comes as one-piece device ready to be filled with a sterile saline solution at the time of surgery.

There are four cuff sizes (3.7cm; 4.0cm; 4.5cm and 5.0cm).

The kit also contains the necessary materials to set up the device before to be implanted. It includes the following disposable items:

- 2 (two) injection non-coring needles: 25G x 5/8”
- 2 (two) 10 ml syringes

ACCESSORIES KIT

The Accessories Kit is delivered to enable the implantation procedure. It includes the following disposable items:

- 1 (one) Measurement tape: to determine the urethral perimeter for correct cuff size selection.
- 30 cm Silicone tubing: protection caps for the instruments used to handle VICTO implant.

ADJUSTMENT KIT

This kit contains the necessary materials to adjust the system internal pressure by adding or extracting sterile saline solution. It includes the following disposable items:

- 2 (two) non-coring injection needles: 25G x 5/8"
- 1 (one) 6 ml syringe

STERILIZATION – SUPPLIED STERILE

All components of the Starter Kit, Accessories Kit, and Adjustment Kit are supplied sterile. Sterility is maintained only if the package is intact and undamaged. Therefore, do not use the product if the package is open or damaged.

SHELF LIFE

For shelf life, refer to the label on the package of individual components or the cardboard box label. Do not use after the expiration date specified in the packaging.

SINGLE USE

All components are intended for SINGLE USE ONLY. DO NOT RE-USE or RE-STERILIZE them, as this can potentially result in compromised device performance and increased risk of inappropriate re-sterilization and cross contamination.

CARE AND STORAGE OF STERILIZED COMPONENTS

In order to protect the packaging integrity and sterile condition of the device and components, it is important that whilst stored, they remain in the cardboard box and are placed on a protected shelf or inside a cabinet. The environment should be clean, dry and near room-temperature.

IMPLANTATION

Prior to implantation of VICTO the surgeon and operating room staff should be familiar with the operating procedure and stages of implantation.

The patient should be fully informed of the procedure and give consent. The surgeon should discuss the possible risks associated with the surgery, anesthesia and device.

PATIENT PREPARATION

Prior to surgery, prophylactic antibiotic therapy should be administered according to the hospital approved procedures.

Once in the operating room, the patient is placed in the lithotomy position, and the abdominal and genital areas are shaved and cleansed with a povidone-iodine solution (provided the patient is not iodine sensitive) or using an approved hospital preoperative scrub procedure.

The patient is prepped and draped for perineal and lower abdominal incision.

EQUIPMENT NEEDED

The implantation of VICTO requires a standard urethroplasty surgical kit as well as some additional supplies. Each surgeon may also have preferences about what equipment he/she requires.

Additional equipment is listed below:

- 12 Fr Foley catheter
- Shod forceps (hemostats)
- Sterile tray
- Sterile 50 ml syringe
- Sterile saline solution
- Umbilical tape

UNPACKING THE COMPONENTS

Prior to unpacking the Starter Kit and Accessories Kit, the integrity of both packages should be inspected – if there is a suspected breach of sterility the device should be returned to the distributor or manufacturer.

Select the Starter Kit after the correct measurement of the urethral size using the measurement tape provided in the Accessories Kit.

The device is supplied with a number of labels detailing the batch number and use by date. These should be retained and placed in the appropriate surgical and patient notes as established in the hospital procedure.

EQUIPMENT PREPARATION

PREPARING HAEMOSTATS

During implantation, it is useful to clamp some of the device tubes. In order to protect the device tubing from damage when clamping, cover the jaws of the hemostat with small sections of silicone tubing provided in Accessories Kit. All the teeth of the hemostat should be covered. The jaws should only be clamped to the first click to prevent excessive pressure on the tubing.

SURGICAL PROCEDURES

BULBAR URETHRAL EXPOSURE

1. Place an indwelling catheter (Foley 12) into the urethra to facilitate identification during the dissection.
2. Make a perineal incision and dissect around the bulbar urethra.
3. Place some umbilical tape around the urethra until the final cuff size is chosen and placed in position.

Figure 3

MEASURE URETHRA SIZE

4. Place the measurement tape around the urethra where the cuff is to be implanted. It should fit snugly without constricting the urethra. NOTE: Before measuring the urethra, take out the indwelling catheter.
5. Select the Starter kit that corresponds to urethra size.

PREPARATION OF THE DEVICE

6. Place the whole device into a sterile tray.
7. Fill the device with 13 ml of sterile saline solution using the non-coring needle 25G x 5/8". (Note: Do not use a different needle than the one provided in the Starter Kit.)
8. Lift the pump upwards, with the filling port at the top, clear of the other components in the tray.
9. Whilst holding the control pump around the filling port between the thumb and index finger, take care when inserting the non-coring needle 25G x 5/8" in-line through the center of the self-sealing port until the tip is just inside the pump chamber.
10. Slowly inject all 13 ml sterile saline solution into the device (perform the process in two steps using the 10ml syringe). The sterile saline solution will enter each of the component parts of the device together with some air. Distribute the sterile saline solution squeezing the top of the control pump.
11. With the control pump raised gently, squeeze the cuff and the PRB. All the air flows into the control pump and vent through the syringe.
12. Finally, remove the non-coring needle 25G x 5/8". Test the system by first pressing the rapid refill valve to re-distribute fluid inside the device, and then squeeze and release the pump twice to check that the cuff deflates. Again, rapidly re-inflate to refill the cuff. A final test is to squeeze and release the pump twice again to deflate the cuff and then check for automatic slow re-inflation. VICTO is now ready for implantation.

ABDOMINAL INCISION

13. Make a lower abdominal incision where the PRB will be placed. Use a blunt dissection to create a space for the balloon.


Figure 4

CUFF PLACEMENT

14. Squeeze the control pump to move the fluid out of the cuff to the PRB and clamp off the tube (cuff – control pump) just below the control pump with a hemostat. Pass the deflated cuff to the perineal site leaving the pump and balloons protruding from the abdominal incision above.

15. The cuff is placed around the urethra and the tab is passed through the loop. Using a hemostat, the button is passed through the hole to secure the cuff in position. The cuff should not be tight (i.e. it should be able to rotate freely around the urethra).

16. The completed cuff insertion should have the appearance shown in the figure, with the connecting tube placed parallel to the urethra by rotating the cuff to the left or right as preferred.

Figure 5

17. Remove the hemostat in order to allow the cuff to refill around the urethra.

INSERTION OF THE PUMP IN THE SCROTUM

18. To place the control pump in the scrotum, use blunt dissection with fingers and push up the scrotal wall on the appropriate side in order that it appears inverted at the abdominal incision. Take the pump and lower it down into a pocket created in the scrotum until it is comfortably situated at the base of the scrotum.

Figure 6

NOTE: The ideal position of the pump in the left or right base of the scrotum should have previously been decided with the patient.

NOTE: If considered necessary, a couple of loose degradable sutures can be placed around the fairing of the tubes leaving the pump to ensure that the pump does not migrate upwards during the initial healing period.

PLACEMENT OF PRB

19. Through the lower abdominal incision place the PRB intraperitoneal.

20. The perineal and abdominal incisions can be closed and dressed.

POST-OPERATIVE CARE

IMMEDIATELY POST-OPERATIVE

An indwelling catheter should remain in place for 24 hours. Antibiotic cover should be

prescribed at the discretion of the surgeon. After 24 hours the indwelling catheter should be removed and the patient should try to void on urge.

If the patient is unable to void without large residual urine, it may be necessary to teach the patient intermittent catheterization in the case of retention.

The patient should continue to use pads or condom drainage to contain leakages.

DISCHARGE

The patient can be discharged 2 to 4 days postoperatively at the discretion of the surgeon and depending on his general health.

Antibiotic cover should be prescribed at the discretion of the surgeon.

An outpatient appointment should be made for 4 to 6 weeks following implantation for pressurization of the device.

PRESSURIZATION PROCEDURE

The VICTO implant is pressurized by injecting sterile saline solution through the self-sealing port in the base of the control pump. The pressure in the device depends upon the amount of fluid injected. This feature allows the device to be set at any desired pressure without having to change the PRB.

EQUIPMENT REQUIRED

The pressurization procedure is best carried out as a sterile procedure in a clean clinical space.

- Sterile saline solution for injection
- VICTO Adjustment Kit
 - 6ml syringe
 - non-coring needles: 25G x 5/8"
- Sterile dressing pack
- Sterile drapes
- Sterile gloves

IMPORTANT

It is vital that the wall of the pump unit is not pierced with the needle. If this occurs, the device may be permanently damaged, necessitating explantation of the whole device.

It is essential that only the non-coring needles 25G x 5/8" provided in the Adjustment Kit be used and that the needle is always aligned parallel to the sides of the pump and penetrates the self-sealing port straight and through the "guard ring." Only the base of the pump with its metal "guard ring" has the ability to be injected safely.

PATIENT PREPARATION

- Explain the procedure to the patient.
- Lay the patient supine and identify the location of the control pump inside the scrotum.
- Cleanse the scrotum in accordance with hospital approved procedure.

NOTE: Local anesthesia around the area may be applied in accordance with hospital approved procedure.

PRESSURIZATION PROCESS

- Fill the syringe with a sterile saline solution for injection.
- Attach the non-coring needle 25G x 5/8" to the syringe.
- With syringe pointed upward, flush 1.0ml of sterile saline solution and expel excess air.
- Fix the self-sealing port of the control pump using two fingers pulling the skin taut over the end of the control pump.
- Push the needle slowly through the skin into the center of the control pump self-sealing port. Just prior to reaching the

control pump, clear the end of the needle by injecting 0.1 ml of sterile saline solution into the tissue space.

- Push the needle into the self-sealing port ensuring that the needle remains parallel with walls of control pump.
- Activate the product injecting 3ml of sterile saline solution.
- Remove the needle slowly.
- Instruct the patient to pump the device about 3 times to redistribute fluid. Rapid re-inflation can be achieved by squeezing the top of the control pump.
- Perform the test of efficacy (See: Post-pressurization continence tests).
- For additional adjustment of the product inject 1ml of sterile saline solution each time.

POST-PRESSURIZATION CONTINENCE TESTS

TESTS OF EFFICACY

These can be done immediately following pressurization process. After showing the patient how to deflate the cuff for voiding, he should be asked to fill his bladder by drinking until a desire to void is reached.

The patient should be asked to give 5 strong coughs, 5 seconds of Valsalva maneuver and if possible 5 squats. The volume leaked can be measured and recorded. The patient should then be asked to void into a flow machine in order that flow rate and volume can be measured and recorded.

LONGER TERM TESTS OF EFFICACY

The patient can be asked to continue wearing pads (or Conveen leg bag) to record pad use and/or leaked volume by weight.

PATIENT EDUCATION

Following pressurization of the device, the surgeon should show the patient how to pump the device, in order to allow urination.

The rapid re-inflate capability can also be demonstrated. If the patient is using intermittent catheterization to empty the bladder, then the importance of pumping before inserting the catheter should be stressed.

A flow rate test should be performed to ensure that the patient is pumping down sufficiently.

An emergency card should be issued to inform healthcare professionals that the patient has an implant and that catheterization should only be performed after pumping down the device to deflate the cuff.

IMAGING

X-RAY VISIBILITY OF VICTO IMPLANT

The VICTO does not contain a radio-opaque solution and therefore, it cannot be imaged using plain radiography (x-ray). However, there are 2 other methods which can be used to image the device in order to investigate function.

MAGNETIC RESONANCE

The device can be imaged using magnetic resonance imaging as described in the paper by Deng et al.

NOTE: see Magnetic Resonance (MR) environment - references.

FLEXIBLE CYSTOMETRY

The function of the cuff can be visualized using flexible cystometry with the cuff inflated and deflated.

ULTRASOUND DIAGNOSTICS

The balloon of VICTO implant can be easily visualized using ultrasound. With appropriate settings, the outline, particularly the highlighted top and bottom of the balloon, can be seen. The diameter of the balloon can be

measured when the cuff is normally inflated. From the diameter of the balloon, and using the following formula, the added volume can be determined. $[(4/3\pi(D/2)^3)-9]$

By using this ultrasound technique, it is also possible to determine the normal functioning of the device during voiding. Pump the sphincter cuff down with about 2-3 good squeezes and notice the enlargement of the PRB. After waiting for about 10-15 minutes, during which time the cuff will be fully re-inflated, the balloon should have returned to approximately the same size.

BIBLIOGRAPHY & REFERENCES:

Knight SL, Susser J, Greenwell T, Mundy AR, Craggs MD. A New Artificial Urinary Sphincter with Conditional Occlusion for Stress Urinary Incontinence: Preliminary Clinical Results. *European Urology*.2006; 50: 574-580.

Deng J, Hall-Craggs MA, Craggs MD, Rickards D, Knight SL, Linney AD, Mundy AR. Three-dimensional MRI of the male urethrae with implanted artificial sphincters: initial results. *The British Journal of Radiology*. 2006;79.

Craggs M, Knight S, Susser J, Kirkham A, Rickards D. Ultrasound imaging as a diagnostic tool for assessing a new artificial urinary sphincter. *Neurourology & Urodynamics*. Proceedings of the International Continence Society, 2006. 174.

García Montes F, Knight SL, Greenwell T, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007; 31(7): 752-758.

Hübner W, Oliver S, Schlarp O. Evaluation of the Flowsecure® Artificial Urinary Sphincter (AUS) Artificial Sphincter: Preliminary results. Proceedings of the International Continence Society, Rotterdam. 2007, 499.

Hübner W. New developments in the area of hydraulic sphincters. *J.Urol.Urogynäkol*. 2007; 14:23-24

García-Montes F, Vicens-Vicens A, Ozonas-Moragues M, Oleza- Simo J, Mora-Salvá A. Surgical implantation of the new Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter in the female bladder neck. *Urol Int*. 2007;79(2):105-10.

García Montes F, Vicens Vicens A, Ozonas Moragues M, Pizá Reus P, Mora Salvá A, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007 Sep;31(8):872-879. [In Spanish].

Khatib-Shahidi, K., Grein M. - Flowsecure® Artificial Urinary Sphincter (AUS) First experiences. Proceedings of NRW Urology, Congress,2008. HELIOS Hospital Schwelm, Germany Poster P60.

Bestard Vallejo JE, García Montes F, Cecchini Rosell L, Mir Maresma C, Morote Robles J. Implantation technique of the artificial urinary sphincter "Flowsecure® Artificial Urinary Sphincter (AUS)" in the bulbar urethra]. *Arch Esp Urol*.2009;62(3):195-200.

Alonso Rodríguez D, Fernández Barranco L, García-Montes F, Mora Salvá A, Ozonas Moragues M. Follow-up of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter. *ActasUrol Esp*. 2009 Oct;33(9):952-5. Review.

García-Montes F. Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter for the treatment of stress urinary incontinence after radical prostatectomy. *Arch Esp Urol*. 2009; 62(10): 845-850. [Spanish].

Fernandez Barranco L, D. Rodriguez A, Garcia Montes F, Mora Salvá A, Ozonas Moragues M. MRI of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter following radical prostatectomy. *Arch. Esp. Urol*.2010; 63 (4): 313-314.

Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. One hundred Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincters. *European Urology Suppl*. 2011; 10(2):309.

Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. Four years experience with the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter: Problems and solution. *Neurourology & Urodynamics*. 2011 Vol Abstract 250.

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BIBLIOGRAPHIE & LITERATURANGABEN:

- Knight SL, Susser J, Greenwell T, Mundy AR, Craggs MD. A New Artificial Urinary Sphincter with Conditional Occlusion for Stress Urinary Incontinence: Preliminary Clinical Results. *European Urology*.2006; 50: 574–580.
- Deng J, Hall-Craggs MA, Craggs MD, Rickards D, Knight SL, Linney AD, Mundy AR. Three-dimensional MRI of the male urethrae with implanted artificial sphincters: initial results. *The British Journal of Radiology*. 2006;79.
- Craggs M, Knight S, Susser J, Kirkham A, Rickards D. Ultrasound imaging as a diagnostic tool for assessing a new artificial urinary sphincter. *Neurourology & Urodynamics*. Proceedings of the International Continence Society, 2006. 174.
- García Montes F, Knight SL, Greenwell T, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007; 31(7): 752-758.
- Hübner W, Oliver S, Schlarp O. Evaluation of the Flowsecure® Artificial Urinary Sphincter (AUS) Artificial Sphincter: Preliminary results. Proceedings of the International Continence Society, Rotterdam. 2007, 499.
- Hübner W. New developments in the area of hydraulic sphincters. *J.Urol.Urogynäkol*. 2007; 14:23-24
- García-Montes F, Vicens-Vicens A, Ozonas-Moragues M, Oleza- Simo J, Mora-Salvá A. Surgical implantation of the new Flowsecu- re® Artificial Urinary Sphincter (AUS) artificial urinary sphincter in the female bladder neck. *Urol Int*. 2007;79(2):105-10.
- García Montes F, Vicens Vicens A, Ozonas Moragues M, Pizá Reus P, Mora Salvá A, Mundy AR, Craggs MD. "Flowsecure® Artificial Urinary Sphincter (AUS)" artificial urinary sphincter: a new adjustable artificial urinary sphincter concept with conditional occlusion for stress urinary incontinence. *Actas Urol Esp*. 2007 Sep;31(8):872-879. [In Spanish].
- Khatib-Shahidi, K., Grein M. - Flowsecure® Artificial Urinary Sphincter (AUS) First experiences. Proceedings of NRW Urology, Congress,2008. HELIOS Hospital Schwelm, Germany Poster P60.
- Bestard Vallejo JE, García Montes F, Cecchini Rosell L, Mir Maresma C, Morote Robles J. Implantation technique of the artificial urinary sphincter "Flowsecure® Artificial Urinary Sphincter (AUS)" in the bulbar urethra]. *Arch Esp Urol*.2009;62(3):195-200.
- Alonso Rodríguez D, Fernández Barranco L, García-Montes F, Mora Salvá A, Ozonas Moragues M. Follow-up of the Flowsecu- re® Artificial Urinary Sphincter (AUS) artificial urinary sphincter. *ActasUrol Esp*. 2009 Oct;33(9):952-5. Review.
- García-Montes F. Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter for the treatment of stress urinary incontinence after radical prostatectomy. *Arch Esp Urol*. 2009; 62(10): 845-850. [Spanish].
- Fernandez Barranco L, D. Rodríguez A, Garcia Montes F, Mora Salvá A, Ozonas Moragues M. MRI of the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter following radical prostatectomy. *Arch. Esp. Urol*.2010; 63 (4): 313-314.
- Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. One hundred Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincters. *European Urology Suppl*. 2011; 10(2):309.
- Alonso Rodríguez D, Fes Ascanio E, Fernández Barranco L, Vicens Vicens A, García-Montes F. Four years experience with the Flowsecure® Artificial Urinary Sphincter (AUS) artificial urinary sphincter: Problems and solution. *Neurourology & Urodynamics*. 2011 Vol Abstract 250.

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Sphinx Medical Limited,